

WIR RX14 Headset Receiver: 2.3 and 2.8 MHz (Stereo)

Instructions For Use & Care

SAFETY INFORMATION

Hearing Safety:



CAUTION!

This product is designed to amplify sounds to a high volume level which could potentially cause hearing damage if used improperly. To protect your hearing and the hearing of others:

1. Make sure the volume is turned down before putting on the earphone or headphone before adjusting the volume to a comfortable level.
2. Set the volume level at the minimum setting that you need to hear.
3. Do not allow children or other unauthorized persons to have access to this product.

BATTERY SAFETY & DISPOSAL



CAUTION!

Do not attempt to recharge disposable batteries, which may explode, release dangerous chemicals, cause burns, or other serious harm to the user or product.

PACEMAKER SAFETY



CAUTION!

1. Before using this product with a pacemaker or other medical device, consult your physician or the manufacturer of your pacemaker or other medical device.
2. If you have a pacemaker or other medical device, make sure that you are using this product in accordance with safety guidelines established by your physician or the pacemaker manufacturer.

RECYCLING INSTRUCTIONS



Help Williams Sound protect the environment! Please take the time to dispose of your equipment properly.

Product Recycling:

Please do NOT dispose of your Williams Sound equipment in the household trash. Please take the equipment to a electronics recycling center; OR, return the product to the factory for proper disposal.

Battery Recycling:

Please do NOT dispose of used batteries in the household trash. Please take the batteries to a retail or community collection point for recycling.

5/10/06


OPERATING INSTRUCTIONS

1. Turn the receiver's Power Switch to ON.
2. Place the headset on your head and adjust the headband size by pulling or pushing on each earpiece.
3. Locate and adjust the individual volume controls for a comfortable listening level.

IMPORTANT

The infrared receiver "eye" is located directly above the battery compartment (see illustration). Do not cover it up or you will lose reception.

Note: To maximize battery life, turn the receiver OFF when it's not in use.



LIMITED WARRANTY

Williams Sound products are engineered, designed and manufactured under carefully controlled conditions to provide you with many years of reliable service.

Williams Sound warrants the Infrared Receiver, RX14 against defects in materials and workmanship for FIVE (5) years. During the first five years from the purchase date, we will promptly repair or replace the Infrared Receiver, RX14.

Earphones, headphones, batteries, cables, carry cases, and all other accessory products carry a 90-day warranty.

WILLIAMS SOUND HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. WILLIAMS SOUND, THEREFORE, DISCLAIMS ALL WARRANTIES NOT SET FORTH ABOVE, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE INFRARED RECEIVER, RX14, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. WILLIAMS SOUND SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL

DAMAGES CAUSED BY ANY USE DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. THE SOLE REMEDY FOR ANY DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT IS REPLACEMENT OF THE PRODUCT. NO PERSON HAS ANY AUTHORITY TO BIND WILLIAMS SOUND TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE INFRARED RECEIVER, RX14. UNAUTHORIZED REPAIRS OR MODIFICATIONS WILL VOID THE WARRANTY.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

If you experience difficulty with your system, call Toll-Free for Customer Assistance:

1-800-843-3544

If it is necessary to return the system for service, your Customer Service Representative will give you a Return Authorization Number (RA) and shipping instructions.

Pack the system carefully and send it to:

Williams Sound Corp.
Attn: Repair Dept.
10321 West 70th Street
Eden Prairie, MN 55344



BATTERY INSTRUCTIONS

This receiver comes standard with disposable AAA Alkaline batteries (BAT 010). Rechargeable AAA NiMH batteries (BAT 022) may also be used. We do not currently offer a charger for this unit. The charging jack is located on the bottom of the left headphone.

EXPECTED BATTERY LIFE

BAT 010 AAA Alkaline batteries	50 hours
BAT 022 AAA NiMH batteries	8 hours per charge

Batteries from other suppliers may provide different operating life.

INSTALLING AND REPLACING BATTERIES

1. The bottom of each battery compartment has a marked “arrow.” Apply pressure to the arrow and pull the compartment open.
2. Remove depleted batteries.
3. To install new batteries, press them into place, *observing proper polarity (+/-)*.
4. Close the battery compartment by inserting the cover top first (tabbed), then snap the bottom into place.
5. When the sound becomes weak or distorted, replace or recharge batteries.
6. If using rechargeable batteries, leave unit charging overnight before operating.

Do not mix batteries and chargers from different suppliers.

**WARNING: DO NOT ATTEMPT TO RECHARGE
NON-RECHARGEABLE BATTERIES!**



WILLIAMS SOUND®

Helping People Hear

10321 West 70th St, Eden Prairie, MN 55344 U.S.A.
800.843.3544 | 952.943.2252 | FAX: 952.943.2174
www.williamsound.com

FRM 441E ©2006 Williams Sound Corp.